

Simultaneous increase in ultrafiltration volume and sodium removal with steady concentration peritoneal dialysis using Carry Life UF

Triomed

Martin Wilkie¹, Jörgen Hegbrant², Charlotte de Leon³, Ola Carlsson³, Olof Heimbürger⁴

1. Sheffield Teaching Hospitals, Sheffield, United Kingdom. 2. Division of Nephrology, Department of Clinical Sciences, Lund University, Lund, Sweden. 3. Triomed AB, Lund, Sweden. 4. Medical Unit Renal Medicine, Karolinska University Hospital and CLINTEC, Karolinska Institute, Stockholm, Sweden.

Introduction

Achieving adequate ultrafiltration (UF) as well as adequate sodium removal is a challenge with peritoneal dialysis (PD), especially during automated PD (APD), as a result of sodium sieving. Carry Life UF is a novel PD system (Figure 1), using steady concentration peritoneal dialysis (SCPD), that is potentially suitable for PD patients requiring enhanced fluid removal. The Carry Life UF device transfers small amounts of dialysate between the patient and the device, while continuously adding glucose to achieve a stable intraperitoneal glucose concentration during the treatment. The aim of this study was to investigate ultrafiltration, sodium removal and glucose efficiency using Carry Life UF.

Methods

- Eight stable PD patients were included in the study
- Subjects were treated with 5-hour Carry Life UF treatments at three different glucose doses (11, 14, 20 g/h)
- An initial fill with 1.5 liters, 1.36% glucose PD solution was used.
- A standard 4-hour 2.27% glucose 2 liters CAPD dwell was used as control
- UF volume, sodium removal, and glucose absorption was calculated as well as glucose efficiency for UF (ml UF volume/g glucose absorbed), and for sodium removal (mmol sodium removed/g glucose absorbed)
- Data expressed as mean \pm SD, statistical analysis using one-way ANOVA, * $p < 0.05$, ** $p < 0.01$, *** $p < 0.001$




Figure 1 Line diagram of the Carry Life UF system: Cyclor, Line set and Glucose bag




Figure 2 Changes in dialysate glucose concentration




Figure 3 Changes in dialysate sodium concentration




Figure 4 Dialytic sodium removal as a function of UF volume; Data from 24 Carry Life UF treatments at three different doses

Treatment	UF volume (mL/dwell)	Sodium removal (mmol/ dwell)	Glucose absorption (g/dwell)	Glucose UF efficiency (mL/g glucose absorbed)	Glucose sodium removal efficiency (mmol/g glucose absorbed)
Control (2.27% CAPD dwell, 4h)	162 \pm 242	21 \pm 33	32 \pm 4.2	5.9 \pm 7.8	0.78 \pm 1.04
Carry Life UF, 5h	11 g/h	646 \pm 256***	86 \pm 27***	43 \pm 9.3**	17.0 \pm 10.6**
	14 g/h	739 \pm 312***	92 \pm 33***	53 \pm 4.4***	14.5 \pm 7.7**
	20 g/h	863 \pm 380***	110 \pm 37***	73 \pm 7.1***	12.4 \pm 6.8*

Table 1 Study results

WCN25-AB-917

Trial registration number:
NCT03724682
(ClinicalTrials.gov)

Reference

Heimbürger O, Hegbrant J, Martus G, Wilkie M, De Leon C, Carlsson O, Johansson AC: Effects of Steady Glucose Concentration Peritoneal Dialysis on Ultrafiltration Volume and Sodium Removal: A Pilot Crossover Trial. Clin J Am Soc Nephrol 2023, 19(2):224-232; 10.2215/CJN.0000000000000342.

The Carry Life UF system is not commercially available

Design of a prospective, multicenter, randomized, cross-over study of the Carry Life UF system in continuous ambulatory peritoneal dialysis (CAPD) patients in the home setting

Martin Wilkie¹, Jörgen Hegbrant², Charlotte de Leon³, Ola Carlsson³, Olof Heimbürger⁴

1. Sheffield Teaching Hospitals, Sheffield, United Kingdom. 2. Division of Nephrology, Department of Clinical Sciences, Lund University, Lund, Sweden. 3. Triomed AB, Lund, Sweden. 4. Medical Unit Renal Medicine, Karolinska University Hospital and CLINTEC, Karolinska Institute, Stockholm, Sweden.

Introduction

Carry Life UF is a novel technology that provides treatment with steady concentration peritoneal dialysis (SCPD) i.e., maintains the glucose concentration in the dialysate for the entire duration of the treatment to obtain an efficient ultrafiltration (UF) and sodium removal. The aim of this study is to assess the efficacy and safety of the Carry Life UF device system (Figure 1) in the home setting. Collection of accurate data in a home environment requires a carefully designed protocol and providing research nursing support for patient training and end-point data collection. (NCT05874804, ClinicalTrials.gov)

Study design

Study type: A prospective, multicenter (Italy, Sweden, UK), randomized, crossover study of adult end-stage kidney disease patients

Study duration: Approximately 12 weeks per participant

Primary endpoint: UF volume, comparing the control CAPD 2.27% glucose dwell with the Carry Life UF treatment, measured at two specific treatments during each arm of the home phase

Secondary endpoints: Adverse events, peritoneal sodium removal, glucose UF efficiency, peak dialysate glucose concentration

Study size: 19 participants required to complete the study

Key inclusion criterion:

A PD prescription of 2–4 CAPD dwells/day unchanged for a minimum of two weeks, with at least one 1.5–2 L, 2.27% glucose day dwell daily

Key exclusion criteria:

- » An episode of peritonitis within the last three months
- » Clinical signs of dehydration
- » Systolic blood pressure < 100 mmHg within the last month

In-clinic phase (3 visits):

- » Peritoneal solute transfer rate classification (PET)
- » Two Carry Life UF treatment (11 g/h or 15 g/h glucose dose with an initial fill of 1.36% glucose PD solution)
- » Hourly dialysate measurements for peak dialysate glucose concentration and dialysate sodium concentration

Home phase:

- » Control arm – 4 weeks on standard PD prescription
- » Carry Life UF arm – 4 weeks with three Carry Life UF treatments per week replacing three 2.27% glucose dwells; the remaining four days, one 2.27% glucose dwell daily is replaced with a 1.36% glucose dwell

Ensuring accurate assessment of study endpoint

Correct determination of UF volumes in PD studies is a challenge. Detailed weighing and sampling instructions are provided to the study centers to ensure accurate data quality. Measures to ensure reliable endpoint data include:

- » Weighing of the CAPD Y-set fill bag performed after the flush and before the fill
- » Weighing and sampling of the CAPD Y-set drain bags performed before the flush-before-fill procedure,
- » Plastic material weight considered in the fluid volume calculations
- » Study centers provided detailed weighing instructions, such as how to place the bags on the scale and weighing the bags with or without overwrap
- » Specially trained study nurses to visit the participants for bag weighing and sampling to obtain reliable data for the endpoint analyses

Economic compensation to participants

Due to the significant burden associated with this device trial and enable a smooth recruitment process, it is important to offer study participants a fair compensation, in accordance with local regulations

Study activities considered when determining compensation:

- » Three in-clinic treatments (whole day visits)
- » Two to five Carry Life UF training sessions
- » Twelve Carry Life UF treatments at home (60 hours plus preparation)
- » 56 days filling out diary (body weight, blood pressure, PD solutions, drain bag weights)
- » Three 24-hour urine collections
- » Two end-of arm visits to the clinic

Conclusion

This study sets out to evaluate a novel PD technology based on SCPD.

The study places a considerable burden on the participants, furthermore, the execution of the study in the home environment entails challenges both in ensuring accurate endpoint data recording and in providing necessary support to the subjects in the use of the technology.

Our protocol has been carefully designed to attend to these issues and to optimize the possibility of the study reaching its stated goals.




Figure 1 Line diagram of the Carry Life UF system




Figure 2 Flow chart presenting the phases of the study.

CLUF = Carry Life UF, PET = peritoneal equilibration test, T = device training